

PRODUCT INFORMATION

Target	IFNAR1
Synonyms	AVP;IFN-alpha-REC;IFNAR;IFNBR;IFRC
Description	Recombinant Human IFNAR1 Protein with C-terminal 6×His tag
Delivery	In Stock
Uniprot ID	P17181
Expression Host	HEK293
Tag	C-6×His Tag
Molecular Characterization	IFNAR1(Lys28-Lys436) 6×His Tag
Molecular Weight	The protein has a predicted molecular mass of 48.0 kDa after removal of the signal peptide. The apparent molecular mass of IFNAR1-His is approximately 70-100 kDa due to glycosylation.
Purity	The purity of the protein is greater than 90% as determined by SDS-PAGE and Coomassie blue staining.
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Storage & Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	The protein encoded by this gene is a type I membrane protein that forms one of the two chains of a receptor for interferons alpha and beta. Binding and activation of the receptor stimulates Janus protein kinases, which in turn phosphorylate several proteins, including STAT1 and STAT2. The protein belongs to the type II cytokine receptor family and functions as an antiviral factor. [provided by RefSeq, Jul 2020]
Usage	Research use only






Figure 1. Human IFNAR1 Protein, His Tag on SDS-PAGE under reducing condition.

Human IFNAR1, His Tagged protein ELISA
0.1 µg of Human IFNAR1, His tagged protein per well




Figure 2. ELISA plate pre-coated by 1 µg/mL (100 µL/well) Human IFNAR1 Protein, His Tag (PME100958) can bind Anti-IFNAR1 Neutralizing antibody BME100117 in a linear range of 0.61-39.06 ng/mL.

