

PRODUCT INFORMATION

Target	SLAMF3
Synonyms	CD229;hly9;mLY9;SLAMF3
Description	Recombinant human SLAMF3 protein with C-terminal 6×His tag
Delivery	In Stock
Uniprot ID	Q9HBG7
Expression Host	HEK293
Tag	C-6×His Tag
Molecular Characterization	SLAMF3(Lys48-Lys454) 6×His tag
Molecular Weight	The protein has a predicted molecular mass of 45.6 kDa after removal of the signal peptide. The apparent molecular mass of SLAMF3-His is approximately 55-70 kDa due to glycosylation.
Purity	The purity of the protein is greater than 90% as determined by SDS-PAGE and Coomassie blue staining.
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Storage&Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Sterility	Products are supplied non-sterile. For cell culture applications, dilute in appropriate medium and sterile-filter (0.22 µm) prior to use.
Background	LY9 belongs to the SLAM family of immunomodulatory receptors (see SLAMF1, MIM 603492) and interacts with the adaptor molecule SAP (SH2D1A, MIM 300490) (Graham et al., 2006 [PubMed 16365421]).
Usage	Research use only
Conjugate	Unconjugated






Figure 1. Human SLAMF3 Protein, His Tag on SDS-PAGE under reducing condition.

DIMABIO CONFIDENTIAL

