

PRODUCT INFORMATION

Target	OX40 Ligand
Synonyms	OX40L;CD252;Tnfsf4
Description	Recombinant mouse OX40 Ligand protein with N-terminal human Fc tag
Delivery	In Stock
Uniprot ID	P43488
Expression Host	HEK293
Tag	N-Human Fc Tag
Molecular Characterization	hFc(Glu99-Ala330) Mouse OX40 Ligand(Gln49-Leu198)
Molecular Weight	The protein has a predicted molecular mass of 43 kDa after removal of the signal peptide. The apparent molecular mass of hFc-mOX40-L is approximately 40-55 kDa due to glycosylation.
Purity	The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Storage&Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Sterility	Products are supplied non-sterile. For cell culture applications, dilute in appropriate medium and sterile-filter (0.22 µm) prior to use.
Background	Cytokine that binds to TNFRSF4. Co-stimulates T-cell proliferation and cytokine production.[UniProtKB/Swiss-Prot Function]
Usage	Research use only
Conjugate	Unconjugated






Figure 1. Mouse OX40 Ligand Protein, hFc Tag on SDS-PAGE under reducing condition.

DIMABIO CONFIDENTIAL

