

PRODUCT INFORMATION

Common Name	LY-3372993;LY 3372993;LY3372993;N3pG 4;N3pG4;N3PG-IV
Conjugate	Unconjugated
Synonyms	AAA;AD1;PN2;ABPP;APPI;CVAP;ABETA;PN-II;preA4;CTFgamma;alpha-sAPP
Applications	N/A
Endotoxin	Less than 1.0 EU/μg by the LAL method. For <1 EU/mg requirements, please contact us for customization.
Recommended Dilutions	N/A
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Host Species	N/A
IgG type	IgG1
Reactivity	Human
Target	APP
Uniprot ID	P05067
Description	Anti-APP(remternetug biosimilar)IgG1) mAb
Delivery	4-6 weeks
Storage&Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Sterility	Products are supplied non-sterile. For cell culture applications, dilute in appropriate medium and sterile-filter (0.22 μm) prior to use.
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

