

PRODUCT INFORMATION

VIS-649; VIS649; VIS 649 **Common Name**

Conjugate Unconjugated

Synonyms APRIL;CD256;TALL2;ZTNF2;TALL-2;TNLG7B;TRDL-1;UNQ383/PRO715

Applications N/A Recommended N/A

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

Host Species IgG type lgG2 Reactivity Human Target **APRIL Uniprot ID** 075888

Description Anti-APRIL(sibeprenlimab biosimilar IgG2) mAb

Delivery 4-6 weeks

Storage & Shipping

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.

Email: info@dimabio.com Website: www.dimabio.com

Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. **Background**

Research use only Usage

