

## **PRODUCT INFORMATION**

**Common Name** BMS-986352;BMS 986352;BMS986352

Conjugate Unconjugated

**Synonyms** BCM;BCMA;CD269;TNFRSF13A

TNFRSF17

**Applications** N/A Recommended N/A

**Dilutions** 

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis Formulation & Reconstitution

for specific instructions of reconstitution.

**Host Species** N/A IgG type lgG1 Reactivity Human **Target BCMA Uniprot ID** Q02223

**Description** Anti-BCMA(ispectamab biosimilar IgG1) mAb

**Delivery** 

Store at -20°C to -80°C for 12 months in

lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store Storage & Shipping

at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient

> Email: info@dimabio.com Website: www.dimabio.com

témperature.

Research grade biosimilar. Not for use in

**Background** therapeutic or diagnostic procedures for humans

or animals.

Usage Research use only

