

PRODUCT INFORMATION

CH 7092230;CH7092230;CH-7092230;RG 6107;RG-6107;RG6107;RO 7092230;RO-7112689;RO7112689;SKY 59;SKY-59;SKY59

Common Name

Conjugate Unconjugated

Synonyms C5D;C5a;C5b;ECLZB;CPAMD4

N/A

Applications N/A Recommended

Dilutions Formulation &

Reconstitution

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis

for specific instructions of reconstitution.

Host Species N/A IgG type lgG1 Reactivity Human **Target** C5

Uniprot ID P01031

Description Anti-C5(crovalimab biosimilar IgG1) mAb

4-6 weeks **Delivery**

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not

intended for use within a month, aliquot and store Storage & Shipping

at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient

Email: info@dimabio.com Website: www.dimabio.com

temperature.

Research grade biosimilar. Not for use in

Background therapeutic or diagnostic procedures for humans

Usage Research use only

