Cat. No. BME100637



PRODUCT INFORMATION

ALXN 1210; ALXN-1210; ALXN1210; ALXN **Common Name** 1810;ALXN-1810;ALXN1810;Ultomiris

Conjugate Unconjugated

Synonyms C5D;C5a;C5b;ECLZB;CPAMD4

Applications N/A Recommended N/A

Dilutions

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

Host Species N/A IgG type IgG2 Reactivity Human C5 **Target Uniprot ID** P01031

Description Anti-C5(ravulizumab biosimilar IgG2) mAb

Delivery 4-6 weeks

> Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not

intended for use within a month, aliquot and store Storage & Shipping at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient

témperature.

Research grade biosimilar. Not for use in

Email: info@dimabio.com Website: www.dimabio.com

Background therapeutic or diagnostic procedures for humans

or animals.

Research use only **Usage**

