Cat. No. BME100727



## **PRODUCT INFORMATION**

**Common Name** CaC-P29;Gohibic;IFX-1;IFX 1;IFX1

**Conjugate** Unconjugated

**Synonyms** C5D;C5a;C5b;ECLZB;CPAMD4

Applications N/A
Recommended N/A
Dilutions

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis

for specific instructions of reconstitution.

Host Species N/A
IgG type IgG4
Reactivity Human
Target C5
Uniprot ID P01031

**Description** Anti-C5(vilobelimab biosimilar□lgG4) mAb

**Delivery** 4-6 weeks

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store

**Storage & Shipping** Intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).

at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient

> Email: info@dimabio.com Website: www.dimabio.com

témperature.

Research grade biosimilar. Not for use in

**Background** therapeutic or diagnostic procedures for humans

or animals.

**Usage** Research use only

