Cat. No. BME100743



PRODUCT INFORMATION

Common Name VTQ-2101;VTQ2101

Conjugate Unconjugated

CRL; GPL; CRL3; GLMR; GLM-R; PLCA2; hGLM-R; IL-31RA; PRO21384; zcytoR17 **Synonyms**

Applications N/A

Recommended Dilutions

N/A

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

Host Species N/A IgG2 IgG type Reactivity Canine

Canine IL31RA **Target Uniprot ID** XP_013963900.1

Anti-Canine IL31RA(nolavetbart biosimilar IgG2) Description

mAb

4-6 weeks **Delivery**

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not

intended for use within a month, aliquot and store Storage & Shipping at -80°C (Avoid repeated freezing and thawing).

Lyophilized proteins are shipped at ambient

Email: info@dimabio.com Website: www.dimabio.com

temperature.

Research grade biosimilar. Not for use in

Background therapeutic or diagnostic procedures for humans

or animals.

Research use only **Usage**

