

PRODUCT INFORMATION

Common Name	CSL324, CSL-324, CSL 324, G-CSFR
Conjugate	Unconjugated
Synonyms	CSF3R; GCSFR; G-CSF receptor; G-CSF-R
Applications	ELISA, Flow Cyt
Endotoxin	Less than 1.0 EU/ μ g by the LAL method. For <1 EU/mg requirements, please contact us for customization.
Recommended Dilutions	ELISA 1:5000-10000, Flow Cyt 1:100
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Host Species	Homo sapiens
IgG type	Human IgG4 - Kappa
Reactivity	Human
Target	CD114
Uniprot ID	Q99062
Description	Anti-CD114(Anumigilimab biosimilar) mAb
Delivery	In Stock
Storage&Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature.
Sterility	Products are supplied non-sterile. For cell culture applications, dilute in appropriate medium and sterile-filter (0.22 μ m) prior to use.
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only



Anti-CD114(Anumigilimab biosimilar) mAb ELISA

0.5 μ g of Human CD114, His tagged protein per well




Figure 1. ELISA plate pre-coated by 5 μ g/mL (100 μ L/well) Human CD114 protein, His Tag (PME100083) can bind Anti-CD114(Anumigilimab biosimilar) mAb (BME100280) in a linear range of 1.37–111.11 ng/mL.




Figure 2. Flow cytometry analysis with 1 μ g/mL Anti-CD114(Anumigilimab biosimilar) mAb (BME100280) on HEK293 cells transfected with Human CD114 protein (Blue histogram) or HEK293 transfected with irrelevant protein (Red histogram).

