Cat. No. BME100423



PRODUCT INFORMATION

Common Name CDX-0159;CDX0159

Conjugate Unconjugated

PBT;SCFR;C-Kit;CD117;MASTC **Synonyms**

Applications N/A Recommended N/A **Dilutions**

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis Formulation & Reconstitution

for specific instructions of reconstitution.

Host Species N/A IgG type lgG1 Reactivity Human **Target** CD117 **Uniprot ID** P10721

Description Anti-CD117(barzolvolimab biosimilar IgG1) mAb

Delivery

Storage & Shipping

Store at -20°C to -80°C for 12 months in

lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).

> Email: info@dimabio.com Website: www.dimabio.com

Lyophilized proteins are shipped at ambient

témperature.

Research grade biosimilar. Not for use in

Background therapeutic or diagnostic procedures for humans

or animals.

Usage Research use only

