

PRODUCT INFORMATION

| | |
|------------------------------|---|
| Common Name | LY3454738;LY-3454738;LY 3454738 |
| Conjugate | Unconjugated |
| Synonyms | OX2R;MOX2R;CD200R;HCRT2 |
| Applications | N/A |
| Recommended Dilutions | N/A |
| Formulation & Reconstitution | Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. |
| Host Species | N/A |
| IgG type | IgG4 |
| Reactivity | Human |
| Target | CD200R1 |
| Uniprot ID | Q8TD46 |
| Description | Anti-CD200R1(ucenprubart biosimilar)IgG4) mAb |
| Delivery | 3-4 weeks |
| Storage & Shipping | Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. |
| Background | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. |
| Usage | Research use only |

