Cat. No. BME100597



PRODUCT INFORMATION

Arzerra; GSK-1841157; GSK1841157; ; GSK **Common Name**

N/A

1841157; Kesimpta; OMB-157; OMB157; OMB 157

Conjugate Unconjugated

Synonyms B1;S7;Bp35;CD20;FMC7;CVID5;LEU-16;MS4A1

Applications N/A

Recommended

Dilutions

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

Host Species N/A IgG type lgG1 Reactivity Human CD20 **Target Uniprot ID** P11836

Description Anti-CD20(ofatumumab biosimilar IgG1) mAb

4-6 weeks **Delivery**

> Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not

intended for use within a month, aliquot and store

Email: info@dimabio.com Website: www.dimabio.com

Storage & Shipping at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient

témperature.

Research grade biosimilar. Not for use in

Background therapeutic or diagnostic procedures for humans

or animals.

Research use only **Usage**

