

PRODUCT INFORMATION

BRIUMV; EMAB-6; EMAB6; EMAB

6;EMAB-603;EMAB603;EMAB 603;LFB-R603;LFB R603;TG-1101;TG1101;RG 1101;TG-20;TG20;TG 20;TGTX-1101;TGTX1101;TGTX 1101 **Common Name**

Conjugate Unconjugated

Synonyms B1;S7;Bp35;CD20;FMC7;CVID5;LEU-16;MS4A1

Applications N/A Recommended N/A

Dilutions

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis Formulation & Reconstitution

for specific instructions of reconstitution.

Host Species N/A IgG type lgG1 Reactivity Human **Target** CD20 **Uniprot ID** P11836

Description Anti-CD20(ublituximab biosimilar IgG1) mAb

4-6 weeks **Delivery**

> Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not

intended for use within a month, aliquot and store Storage & Shipping

at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient

Email: info@dimabio.com Website: www.dimabio.com

temperature.

Research grade biosimilar. Not for use in

Background therapeutic or diagnostic procedures for humans

Usage Research use only

