

PRODUCT INFORMATION

Common Name	MOR 202;MOR202;MOR-202;MOR-03087;MOR 03087;MOR03087;MOR 3087;MOR-3087;MOR3087;TJ-202;TJ202;TJ 202
Conjugate	Unconjugated
Synonyms	CD38;ADPRC1;cADPR1;ADPRC 1;T10;cADPr hydrolase 1
Applications	N/A
Endotoxin	Less than 1.0 EU/µg by the LAL method. For <1 EU/mg requirements, please contact us for customization.
Recommended Dilutions	N/A
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Host Species	N/A
IgG type	IgG1
Reactivity	Human
Target	CD38
Uniprot ID	P28907
Description	Anti-CD38(felzartamab biosimilar) IgG1) mAb
Delivery	4-6 weeks
Storage&Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Sterility	Products are supplied non-sterile. For cell culture applications, dilute in appropriate medium and sterile-filter (0.22 µm) prior to use.
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

