

PRODUCT INFORMATION

| | |
|---|---|
| Common Name | MOR 202;MOR202;MOR-202;MOR-03087;MOR 03087;MOR03087;MOR 3087;MOR-3087;MOR3087;TJ-202;TJ202;TJ 202 |
| Conjugate | Unconjugated |
| Synonyms | CD38;ADPRC1;cADPR1;ADPRC 1;T10;cADPr hydrolase 1 |
| Applications | N/A |
| Endotoxin | Less than 1.0 EU/μg by the LAL method. For <1 EU/mg requirements, please contact us for customization. |
| Recommended Dilutions | N/A |
| Formulation & Reconstitution | Lyophilized from sterile PBS, pH 7.4. Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. |
| Host Species | N/A |
| IgG type | IgG1 |
| Reactivity | Human |
| Target | CD38 |
| Uniprot ID | P28907 |
| Description | Anti-CD38(felzartamab biosimilar) IgG1) mAb |
| Delivery | 4-6 weeks |
| Storage&Shipping | Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. |
| Background | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. |
| Usage | Research use only |

