

PRODUCT INFORMATION

Common Name TBC

Conjugate Unconjugated

CD39;SPG64;ATP-DPH;ATPDase;NTPDase-1;ENTPD1 **Synonyms**

Applications N/A

Recommended **Dilutions**

Background

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for Formulation & Reconstitution

specific instructions of reconstitution.

Host Species N/A IgG type IgG4 Reactivity Human **CD39 Target Uniprot ID** P49961

Description Anti-CD39(eltivutabart biosimilar IgG4) mAb

4-6 weeks **Delivery**

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid Storage & Shipping

repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.

Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or

> Email: info@dimabio.com Website: www.dimabio.com

animals.

Usage Research use only

