

PRODUCT INFORMATION

ChAglyCD3;GSK-2136525;GSK2136525;GSK 2136525;TRX4;TRX-4;TRX 4 **Common Name**

Conjugate Unconjugated

Synonyms T3E;TCRE;IMD18;CD3epsilon

Applications N/A

Recommended

N/A Dilutions

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

Host Species N/A IgG type IgG1 Reactivity Human CD3E **Target Uniprot ID** P07766

Description Anti-CD3E(otelixizumab biosimilar IgG1) mAb

Delivery 4-6 weeks

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not

intended for use within a month, aliquot and store Storage & Shipping at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient

témperature.

Research grade biosimilar. Not for use in

Email: info@dimabio.com Website: www.dimabio.com

Background therapeutic or diagnostic procedures for humans

or animals.

Research use only **Usage**

