Cat. No. BME100420



PRODUCT INFORMATION

INCA-34176;INCA34176;SNDX-6352;SNDX6352;UCB 6352;UCB6352;UCB-6352 Common Name

Coniugate Unconjugated

 ${\sf FMS;CSFR;FIM2;GPSC;HDLS;C-FMS;CD115;HDLS1;CSF-1R;BANDDOS;M-CSF-R}$ **Synonyms**

Applications Recommended N/A **Dilutions**

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

Host Species IgG type IgG4 Reactivity Human Target CSF1R **Uniprot ID** P07333

Description Anti-CSF1R(axatilimab biosimilar | IgG4) mAb

Delivery

Storage & Shipping

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.

Email: info@dimabio.com Website: www.dimabio.com

Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. **Background**

Research use only Usage



