

PRODUCT INFORMATION

Common Name	INCA-34176;INCA34176;SNDX-6352;SNDX6352;UCB 6352;UCB6352;UCB-6352
Conjugate	Unconjugated
Synonyms	FMS;CSFR;FIM2;GPSC;HDLS;C-FMS;CD115;HDLS1;CSF-1R;BANDDOS;M-CSF-R
Applications	N/A
Endotoxin	Less than 1.0 EU/ μ g by the LAL method. For <1 EU/mg requirements, please contact us for customization.
Recommended Dilutions	N/A
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Host Species	N/A
IgG type	IgG4
Reactivity	Human
Target	CSF1R
Uniprot ID	P07333
Description	Anti-CSF1R(axatilimab biosimilar \square IgG4) mAb
Delivery	4-6 weeks
Storage&Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Sterility	Products are supplied non-sterile. For cell culture applications, dilute in appropriate medium and sterile-filter (0.22 μ m) prior to use.
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

