Cat. No. BME100707



## **PRODUCT INFORMATION**

CP-675;CP675;CP 675;CP-675206;CP 675206;CP675206;IMJUDO;Ticilimumab **Common Name** 

Conjugate Unconjugated

**Synonyms** CD;GSE;GRD4;ALPS5;CD152;CTLA-4;IDDM12;CELIAC3

**Applications** N/A Recommended N/A Dilutions

Formulation & Reconstitution Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

**Host Species** IgG type lgG2 Reactivity Human **Target** CTLA-4 **Uniprot ID** P16410

Description Anti-CTLA-4(tremelimumab biosimilar IgG2) mAb

**Delivery** 4-6 weeks

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins

Email: info@dimabio.com Website: www.dimabio.com

Storage & Shipping

are shipped at ambient temperature.

Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. **Background** 

Usage Research use only

