

PRODUCT INFORMATION

|                              |   |
|------------------------------|---|
| Common Name                  | AGEN-1884;AGEN1884;AGEN 1884;RebmAb-600;RebmAb 600;RebmAb600;UGN-301;UGN301;UGB 301   |
| Conjugate                    | Unconjugated  |
| Synonyms                     | CD;GSE;GRD4;ALPS5;CD152;CTLA-4;IDDM12;CELIAC3   |
| Applications                 | N/A   |
| Recommended Dilutions        | N/A   |
| Formulation & Reconstitution | Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.  |
| Host Species                 | N/A   |
| IgG type                     | IgG1  |
| Reactivity                   | Human   |
| Target                       | CTLA-4  |
| Uniprot ID                   | P16410  |
| Description                  | Anti-CTLA-4(zalifrelimab biosimilar IgG1) mAb   |
| Delivery                     | 3-4 weeks   |
| Storage & Shipping           | Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. |
| Background                   | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.   |
| Usage                        | Research use only   |

