

PRODUCT INFORMATION

Common Name GMA-301;GMA301;GMA 301;Getagozumab

Conjugate Unconjugated

ETA;ET-A;ETAR;ETRA;MFDA;ETA-R;hET-AR **Synonyms**

Applications N/A Recommended N/A

Dilutions

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis Formulation & Reconstitution

for specific instructions of reconstitution.

Host Species N/A IgG type IgG4 Reactivity Human **Target EDNRA Uniprot ID** P25101

Description Anti-EDNRA(patecibart biosimilar IgG4) mAb

Delivery 4-6 weeks

Store at -20°C to -80°C for 12 months in

lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store Storage & Shipping

at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient

témperature.

Research grade biosimilar. Not for use in

Background therapeutic or diagnostic procedures for humans

or animals.

Usage Research use only



Email: info@dimabio.com Website: www.dimabio.com

