

PRODUCT INFORMATION

Common Name JMT101;JMT-101 Conjugate Unconjugated

Synonyms ERBB;ERRP;HER1;mENA;ERBB1;NNCIS;PIG61;NISBD2

Applications N/A Recommended N/A **Dilutions**

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation &

Reconstitution

Host Species N/A IgG type lgG1 Reactivity Human

Target EGFR Uniprot ID P00533

Description Anti-EGFR(becotatug biosimilar IgG1) mAb

Delivery 4-6 weeks

Store at -20 $^{\circ}$ C to -80 $^{\circ}$ C for 12 months in lyophilized form. After reconstitution, if not intended for use

within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. Storage & Shipping

Research grade biosimilar. Not for use in therapeutic **Background**

or diagnostic procedures for humans or animals.

Email: info@dimabio.com Website: www.dimabio.com

Usage Research use only

