

## **PRODUCT INFORMATION**

DS-992;DS992;S-95026;S95206;Sym-004 **Common Name** (Combination with Modotuximab); Sym004

(Combination with Modotuximab)

Conjugate Unconjugated

ERBB;ERRP;HER1;mENA;ERBB1;NNCIS;PIG61;NISBD2 **Synonyms** 

**Applications** N/A Recommended

N/A Dilutions

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for Formulation & Reconstitution specific instructions of reconstitution.

**Host Species** N/A IgG type lgG1 Reactivity Human **Target FGFR Uniprot ID** P00533

Description Anti-EGFR(futuximab biosimilar□lgG1) mAb

**Delivery** 4-6 weeks

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins Storage & Shipping

are shipped at ambient temperature.

Research grade biosimilar. Not for use in therapeutic **Background** 

or diagnostic procedures for humans or animals.

Email: info@dimabio.com Website: www.dimabio.com

Usage Research use only

