

## **PRODUCT INFORMATION**

**Common Name** ABX-EGF;E7.6.3;Panitumab;Vectibix

Conjugate Unconjugated

**Synonyms** ERBB; ERRP; HER1; mENA; ERBB1; NNCIS; PIG61; NISBD2

**Applications** N/A

Recommended **Dilutions** 

N/A

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

**Host Species** N/A IgG type IgG2 Reactivity Human **Target EGFR Uniprot ID** P00533

**Description** Anti-EGFR(panitumumab biosimilar IgG2) mAb

**Delivery** 4-6 weeks

Store at -20  $^{\circ}$ C to -80  $^{\circ}$ C for 12 months in lyophilized form. After reconstitution, if not intended for use

within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. Storage & Shipping

Research grade biosimilar. Not for use in therapeutic **Background** 

or diagnostic procedures for humans or animals.

Email: info@dimabio.com Website: www.dimabio.com

**Usage** Research use only

