

PRODUCT INFORMATION

AB-023;AB 023;AB023;AB022;AB **Common Name**

N/A

022;AB-022;Xisomab 3G3

Conjugate Unconjugated

Synonyms FXI;PTA

Applications N/A

Recommended

Dilutions

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

Host Species N/A IgG type IgG4 Reactivity Human F11 **Target Uniprot ID** P03951

Description Anti-F11(gruticibart biosimilar ☐ IgG4) mAb

Delivery 4-6 weeks

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not

intended for use within a month, aliquot and store Storage & Shipping at -80°C (Avoid repeated freezing and thawing).

Lyophilized proteins are shipped at ambient

Email: info@dimabio.com Website: www.dimabio.com

témperature.

Research grade biosimilar. Not for use in

Background therapeutic or diagnostic procedures for humans

or animals.

Research use only **Usage**

