

PRODUCT INFORMATION

Common Name PYX-201;PYX201;PYX 201

Conjugate Unconjugated

FN;CIG;FNZ;MSF;ED-B;FINC;GFND;LETS;GFND2;SMDCF **Synonyms**

Applications N/A Recommended

N/A **Dilutions**

Lyophilized from sterile PBS, pH 7.4. Normally 5 % -8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for Formulation & Reconstitution

specific instructions of reconstitution.

Host Species N/A IgG type lgG1 Reactivity Human FN1 **Target** P02751 **Uniprot ID**

Description Anti-FN1(micvotabart biosimilar IgG1) mAb

Delivery 4-6 weeks

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins

Storage & Shipping

are shipped at ambient temperature.

Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. **Background**

Email: info@dimabio.com Website: www.dimabio.com

Usage Research use only

