

PRODUCT INFORMATION

Common Name AV-380;AV 380;AV380

Conjugate Unconjugated

HG;PDF;MIC1;PLAB;MIC-1;NAG-1;PTGFB;GDF-15 **Synonyms**

Applications N/A Recommended N/A

Dilutions

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis Formulation & Reconstitution

for specific instructions of reconstitution.

Host Species N/A IgG type lgG1 Reactivity Human **Target** GDF15 **Uniprot ID** Q99988

Description Anti-GDF15(rilogrotug biosimilar IgG1) mAb

Delivery 4-6 weeks

Store at -20°C to -80°C for 12 months in

lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Storage & Shipping

Lyophilized proteins are shipped at ambient

Email: info@dimabio.com Website: www.dimabio.com

témperature.

Research grade biosimilar. Not for use in

Background therapeutic or diagnostic procedures for humans

or animals.

Usage Research use only

