

PRODUCT INFORMATION

Clone ID	OTI2A4 (formerly 2A4)
Target	GST tag
Synonyms	ANTI-GST; ANTIBODY FOR GST; ANTI GST; GST AB; GST ANTIBODY
Host Species	Mouse
Description	Mouse monoclonal eGFP antibody, clone OTI5A2
Delivery	1-2week
Uniprot ID	N/A
IgG type	IgG1
Clonality	Monoclonal
Reactivity	N/A
Applications	WB
Recommended Dilutions	WB: 1:4000
Purification	N/A
Endotoxin	Less than 1.0 EU/ μ g by the LAL method. For <1 EU/mg requirements, please contact us for customization.
Formulation & Reconstitution	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Storage&Shipping	Store at -20°C.
Sterility	Products are supplied non-sterile. For cell culture applications, dilute in appropriate medium and sterile-filter (0.22 μ m) prior to use.
Background	N/A
Usage	Research use only
Conjugate	Unconjugated






Figure 1. Anti-GST immunoblot analysis of six different HEK293T over-expression lysates containing different GST tagged recombinant proteins, and one GST-DDK bacterial lysate produced in E.coli. at 1:5000 dilution

DIMABIO CONFIDENTIAL

