

PRODUCT INFORMATION

Clone ID	N/A
Target	HEL Rabbit IgG1
Synonyms	N/A
Host Species	Rabbit
Description	Anti-HEL Rabbit IgG1 Isotype Control antibody
Delivery	4-6 weeks
Uniprot ID	N/A
IgG type	N/A
Clonality	Monoclonal
Reactivity	Rabbit
Applications	N/A
Recommended Dilutions	N/A
Purification	Purified from cell culture supernatant by affinity chromatography
Endotoxin	Less than 1.0 EU/μg by the LAL method. For <1 EU/mg requirements, please contact us for customization.
Formulation & Reconstitution	Powder
Storage & Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).
Sterility	Products are supplied non-sterile. For cell culture applications, dilute in appropriate medium and sterile-filter (0.22 μm) prior to use.
Background	N/A
Usage	Research use only
Conjugate	Unconjugated

DIMABIO CONFIDENTIAL

