

## **PRODUCT INFORMATION**

**Common Name** XMT-2056; XMT2056; XMT 2056

N/A

Conjugate Unconjugated

NEU;NGL;HER2;TKR1;CD340;HER-2;VSCN2;MLN 19;MLN-19;c-ERB2;c-ERB-2;HER-2/neu;p185(erbB2) **Synonyms** 

**Applications** N/A

Recommended **Dilutions** 

Lyophilized from sterile PBS, pH 7.4. Normally 5 %

- 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for Formulation & Reconstitution

specific instructions of reconstitution.

**Host Species** N/A lgG1 IgG type Reactivity Human **Target** Her2

P04626 **Uniprot ID** 

Anti-Her2(calotatug biosimilar IgG1) mAb **Description** 

**Delivery** 4-6 weeks

Store at -20°C to -80°C for 12 months in

lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Storage & Shipping

Lyophilized proteins are shipped at ambient

Email: info@dimabio.com Website: www.dimabio.com

temperature.

**Background** 

Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans

**Usage** Research use only

