

PRODUCT INFORMATION

| | |
|---|---|
| Common Name | R-1273;R1273;R 1273;RG-1273;RG 1273;RG1273;RO-4368451;RO4368451;RO 4368451 |
| Conjugate | Unconjugated |
| Synonyms | NEU;NGL;HER2;TKR1;CD340;HER-2;VSCN2;MLN 19;MLN-19;c-ERB2;c-ERB-2;HER-2/neu;p185(erbB2) |
| Applications | N/A |
| Endotoxin | Less than 1.0 EU/μg by the LAL method. For <1 EU/mg requirements, please contact us for customization. |
| Recommended Dilutions | N/A |
| Formulation & Reconstitution | Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. |
| Host Species | N/A |
| IgG type | IgG1 |
| Reactivity | Human |
| Target | Her2 |
| Uniprot ID | P04626 |
| Description | Anti-Her2(pertuzumab biosimilar) IgG1) mAb |
| Delivery | 4-6 weeks |
| Storage&Shipping | Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. |
| Sterility | Products are supplied non-sterile. For cell culture applications, dilute in appropriate medium and sterile-filter (0.22 μm) prior to use. |
| Background | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. |
| Usage | Research use only |

