Cat. No. BME100715



## **PRODUCT INFORMATION**

YYB101;YYB-101;YYB 101;CLM-101;CLM 101;CLM101;NOV 1105;NOV-1105;NOV1105 **Common Name** 

Conjugate Unconjugated

SF;HGFB;HPTA;F-TCF;DFNB39 **Synonyms** 

**Applications** N/A Recommended N/A Dilutions

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

**Host Species** N/A IgG type IgG4 Reactivity Human **HGF Target Uniprot ID** P14210

**Description** Anti-HGF(umikibart biosimilar IgG4) mAb

4-6 weeks **Delivery** 

> Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not

intended for use within a month, aliquot and store Storage & Shipping at -80°C (Avoid repeated freezing and thawing).

Lyophilized proteins are shipped at ambient

Email: info@dimabio.com Website: www.dimabio.com

témperature.

Research grade biosimilar. Not for use in

**Background** therapeutic or diagnostic procedures for humans

or animals.

Research use only **Usage** 

