

PRODUCT INFORMATION

ABT-308;ABT308;CC 93538;CC93538;CC-93538;RPC 4046;RPC4046;RPC-4046 **Common Name**

Conjugate Unconjugated **Synonyms** P600;IL-13;NC30

Applications N/A

Recommended **Dilutions**

N/A

Formulation & Reconstitution Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Host Species N/A IgG type IgG1 Reactivity Human **Target** IL13

P35225 **Uniprot ID Description** Anti-IL13(cendakimab biosimilar IgG1) mAb

Delivery 4-6 weeks

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not

intended for use within a month, aliquot and store Storage & Shipping at -80°C (Avoid repeated freezing and thawing).

Lyophilized proteins are shipped at ambient


temperature.

Research grade biosimilar. Not for use in **Background**

therapeutic or diagnostic procedures for humans

or animals.

Usage Research use only



Email: info@dimabio.com Website: www.dimabio.com

