

PRODUCT INFORMATION

Common Name	AIN457
Synonyms	CTLA-8;CTLA8;IL-17;IL-17A;IL17
Applications	ELISA, Flow Cyt
Recommended Dilutions	ELISA 1:5000-10000, Flow Cyt 1:100
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Host Species	Homo sapiens
IgG type	IgG1
Reactivity	Human
Target	IL17A
Uniprot ID	Q16552
Description	Anti-IL17A(secukinumab biosimilar) mAb
Delivery	In Stock
Storage & Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only






Figure 1. Flow cytometry analysis under cell membrane permeable condition with 1 μ g/mL Anti-IL17A (secukinumab biosimilar) mAb (BME100137) on Expi293 cells transfected with Human IL17A protein (Blue histogram) or Expi293 transfected with irrelevant protein (Red histogram).

Anti-IL17A(secukinumab biosimilar) mAb ELISA
0.2 μ g of Human IL17A, hFc tagged protein per well




Figure 2. ELISA plate pre-coated by 2 μ g/mL (100 μ L/well) Human IL17A Protein, hFc Tag (PME100884) can bind Anti-IL17A(secukinumab biosimilar) mAb (BME100137) in a linear range of 6.40-160 ng/mL. In order to specifically detect BME100137, mouse anti-human Fab-specific antibody was used as detection antibody.

