

PRODUCT INFORMATION

Common Name	AMG827, KHK4827
Conjugate	Unconjugated
Synonyms	CD217;CDw217;IL-17RA;IL17R;CANDF5;hIL-17R
Applications	ELISA; Flow Cyt
Recommended Dilutions	ELISA 1:5000-10000; Flow Cyt 1:100
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Host Species	Homo sapiens
IgG type	Human IgG2 - Kappa
Reactivity	Human
Target	IL17RA
Uniprot ID	Q96F46
Description	Anti-IL17RA (brodalumab biosimilar) mAb
Delivery	In Stock
Storage & Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
DIMA Disclaimer	All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement.



Anti-IL17RA (brodalumab biosimilar) mAb ELISA

0.2 µg of Human IL17RA, His Tagged protein per well




Figure 1. ELISA plate pre-coated by 2 µg/ml (100 µl/well) Human IL17RA, His tagged protein ([getskuurl sku="PME100468"]) can bind Anti-IL17RA Neutralizing antibody (BME100038) in a linear range of 3.2-80 ng/ml.




Figure 2. HEK293 cell line transfected with irrelevant protein (B) and human IL17RA (A) were surface stained with Anti-IL17RA mAb 1 µg/ml (brodalumab) followed by Alexa 488-conjugated anti-human IgG secondary antibody.

