

PRODUCT INFORMATION

Common Name	CERC-007, AEVI-007, AVTX-007, CERC007
Synonyms	IGIF; IL1F4;IL-18
Conjugate	Unconjugated
Applications	N/A
Recommended Dilutions	N/A
Endotoxin	Less than 1.0 EU/ μ g by the LAL method. For <1 EU/mg requirements, please contact us for customization.
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Host Species	Homo sapiens
IgG type	Human IgG1(K214R,L234F,L235E,P331S) - Kappa
Reactivity	Human
Target	IL18
Uniprot ID	Q14116
Description	Anti-IL18(camoteskimab biosimilar) mAb
Delivery	In Stock
Storage&Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature.
Sterility	Products are supplied non-sterile. For cell culture applications, dilute in appropriate medium and sterile-filter (0.22 μ m) prior to use.
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only






Figure 1. Anti-IL18(camoteskimab biosimilar) mAb on SDS-PAGE under reducing (R) and non-reducing (NR) conditions.

DIMABIO CONFIDENTIAL

