

PRODUCT INFORMATION

Common Name	S-78989;S78989;S 78989;VPM 087;VPM087;VPM-087;VPM87;VPM-87;XMA-005.2;XMA-0052;XMA52;XOMA-052;XOMA052;XOMA52;XOMA-52
Conjugate	Unconjugated
Synonyms	IL-1;IL1F2;IL1beta;IL1-BETA
Applications	N/A
Recommended Dilutions	N/A
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Host Species	N/A
IgG type	IgG2
Reactivity	Human
Target	IL1B
Uniprot ID	P01584
Description	Anti-IL1B(gevokizumab biosimilar[IgG2) mAb
Delivery	4-6 weeks
Storage&Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

DIMABIO CONFIDENTIAL

