

PRODUCT INFORMATION

Common Name LY 3074828; LY-3074828; LY3074828

Conjugate Unconjugated

P19;SGRF;IL-23;IL-23A;IL23P19 **Synonyms**

Applications N/A Recommended N/A **Dilutions**

Formulation & Reconstitution Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis

for specific instructions of reconstitution.

Host Species N/A IgG type IgG4 Reactivity Human **Target** IL23A **Uniprot ID** Q9NPF7

Description Anti-IL23A(mirikizumab biosimilar IgG4) mAb

Delivery 4-6 weeks

Store at -20°C to -80°C for 12 months in

lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store Storage & Shipping

at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient

témperature.

Research grade biosimilar. Not for use in

Background therapeutic or diagnostic procedures for humans

or animals.

Usage Research use only







