Cat. No. BME100653



PRODUCT INFORMATION

ABBV 066;ABBV066;ABBV 066;ABBV **Common Name**

66;ABBV-66;ABBV66;BI-655066;BI655066;BI

655066;SKYRIZ

Conjugate Unconjugated

Synonyms P19;SGRF;IL-23;IL-23A;IL23P19

Applications N/A

Recommended

Dilutions

N/A

Formulation & Reconstitution Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Host Species N/A IgG type lgG1

Reactivity Human **Target** IL23A Q9NPF7

Uniprot ID

Description Anti-IL23A(risankizumab biosimilar IgG1) mAb

Delivery 4-6 weeks

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not

intended for use within a month, aliquot and store Storage & Shipping at -80°C (Avoid repeated freezing and thawing).

Lyophilized proteins are shipped at ambient

temperature.

Research grade biosimilar. Not for use in **Background**

therapeutic or diagnostic procedures for humans

Email: info@dimabio.com Website: www.dimabio.com

or animals.

Usage Research use only



