

PRODUCT INFORMATION

Common Name RG 6292;RG-6292;RG6962;RO 7296682;RO-7296682;RO7296682

Conjugate Unconjugated

Synonyms p55;CD25;IL2R;IMD41;TCGFR;IDDM10;IL2-RA;IL-2-RA

Applications N/A Recommended

Dilutions N/A

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for

specific instructions of reconstitution.

Host Species N/A
IgG type IgG1
Reactivity Human
Target IL2RA
Uniprot ID P01589

Description Anti-IL2RA(vopikitug biosimilar□lgG1) mAb

Delivery 4-6 weeks

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid within a month, aliquot than store at -80°C (Avoid within a month, aliquot than store).

repeated freezing and thawing). Lyophilized proteins

Email: info@dimabio.com Website: www.dimabio.com

are shipped at ambient temperature. Research grade biosimilar. Not for use in

Background therapeutic or diagnostic procedures for humans or

animals.

Usage Research use only

