

PRODUCT INFORMATION

LY-3375880;LY3375880;LY **Common Name** 3375880;ZB-880;ZB880;ZB 880

Conjugate Unconjugated

Synonyms DVS27;IL1F11;NF-HEV;NFEHEV;C9orf26

Applications N/A Recommended N/A

Dilutions

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

Host Species N/A IgG type IgG4 Reactivity Human IL33 **Target Uniprot ID** 095760

Description Anti-IL33(torudokimab biosimilar IgG4) mAb

4-6 weeks **Delivery**

> Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not

intended for use within a month, aliquot and store Storage & Shipping at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient

témperature.

Research grade biosimilar. Not for use in

Email: info@dimabio.com Website: www.dimabio.com

Background therapeutic or diagnostic procedures for humans

or animals.

Research use only **Usage**

