Cat. No. BME100634



PRODUCT INFORMATION

Common Name CBP 201;CBP-201; CBP201

Conjugate Unconjugated

CD124;IL4RA;IL-4RA **Synonyms**

Applications N/A Recommended N/A **Dilutions**

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis Formulation & Reconstitution

for specific instructions of reconstitution.

Host Species N/A IgG type IgG4 Reactivity Human **Target** IL4RA **Uniprot ID** P24394

Description Anti-IL4RA(rademikibart biosimilar IgG4) mAb

Delivery 3-4 weeks

Store at -20°C to -80°C for 12 months in

lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store Storage & Shipping at -80°C (Avoid repeated freezing and thawing).

Lyophilized proteins are shipped at ambient

Email: info@dimabio.com Website: www.dimabio.com

témperature.

Research grade biosimilar. Not for use in

Background therapeutic or diagnostic procedures for humans

or animals.

Usage Research use only

