

PRODUCT INFORMATION

Common Name	CDP-835;CDP835;CDP 835;CEP-38072;CEP38072;CEP 38072;CINQAERO;CINQAIR;Cinquil;CTx-55700;CTx55700;Ctx 55700;DCP-835;DCP835;DCP 835;SCH 55700;SCH55600;SCH-55700;TRFK-5;TRFK 5;TRFK5
Conjugate	Unconjugated
Synonyms	EDF;TRF;IL-5;Interleukin-5
Applications	N/A
Endotoxin	Less than 1.0 EU/µg by the LAL method. For <1 EU/mg requirements, please contact us for customization.
Recommended Dilutions	N/A
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Host Species	N/A
IgG type	IgG4
Reactivity	Human
Target	IL5
Uniprot ID	P05113
Description	Anti-IL5(reslizumab biosimilar) IgG4) mAb
Delivery	4-6 weeks
Storage&Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Sterility	Products are supplied non-sterile. For cell culture applications, dilute in appropriate medium and sterile-filter (0.22 µm) prior to use.
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

