

PRODUCT INFORMATION

PF-04236921;PF-4236921;PF 04236921;PF 4236921;PF04236921;PF4236921;TOUR-006;TOUR006;TOUR **Common Name**

Conjugate Unconjugated

CDF;HGF;HSF;BSF2;IL-6;BSF-2;IFNB2;IFN-beta-2 **Synonyms**

Applications Recommended N/A **Dilutions**

Formulation & Reconstitution Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of

reconstitution.

Host Species N/A IgG2 IgG type Reactivity Human IL6 **Target Uniprot ID** P05231

Description Anti-IL6(pacibekitug biosimilar ☐IgG2) mAb

Delivery

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient

> Email: info@dimabio.com Website: www.dimabio.com

Storage & Shipping

temperature.

Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. **Background**

Research use only Usage

