

PRODUCT INFORMATION

RZ-358;X 358;XOMA **Common Name**

247;XOMA-247;XOMA247;XOMA

358;XOMA-358;XOMA358;XPA.15.247

Conjugate Unconjugated

Synonyms HHF5;CD220

Applications N/A

Recommended

N/A **Dilutions**

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis Formulation & Reconstitution

for specific instructions of reconstitution.

Host Species N/A IgG type IgG2 Reactivity Human **Target INSR**

P06213 **Uniprot ID**

Anti-INSR(ersodetug biosimilar IgG2) mAb **Description**

Delivery 4-6 weeks

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not

intended for use within a month, aliquot and store Storage & Shipping

at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient

temperature.

Research grade biosimilar. Not for use in **Background**

therapeutic or diagnostic procedures for humans

Email: info@dimabio.com Website: www.dimabio.com

or animals.

Usage Research use only

