

## **PRODUCT INFORMATION**

**Common Name TBC** 

Conjugate Unconjugated

kidney associated antigen 1;kidney-associated **Synonyms** antigen 1;MGC78738;RU2 antisense gene

protein; RÚ2ÁS

**Applications** N/A

Recommended

**Dilutions** 

N/A

Formulation & Reconstitution Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

**Host Species** N/A IgG type IgG1 Reactivity Human **Target** KAAG1 Q9UBP8 **Uniprot ID** 

Anti-KAAG1(olintatug biosimilar IgG1) mAb **Description** 

**Delivery** 4-6 weeks

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not

intended for use within a month, aliquot and store Storage & Shipping at -80°C (Avoid repeated freezing and thawing).

Lyophilized proteins are shipped at ambient

temperature.

Research grade biosimilar. Not for use in **Background** 

therapeutic or diagnostic procedures for humans

Email: info@dimabio.com Website: www.dimabio.com

or animals.

**Usage** Research use only

