



PRODUCT INFORMATION

|                              |                                                                                                                                                                                                                                                                        |
|------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Common Name                  | LAG-3,CD223                                                                                                                                                                                                                                                            |
| Conjugate                    | Unconjugated                                                                                                                                                                                                                                                           |
| Synonyms                     | BMS-986016;ONO-4482                                                                                                                                                                                                                                                    |
| Applications                 | ELISA; Flow Cyt                                                                                                                                                                                                                                                        |
| Recommended Dilutions        | ELISA 1:5000-10000; Flow Cyt 1:100                                                                                                                                                                                                                                     |
| Formulation & Reconstitution | Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.                                                                       |
| Host Species                 | Homo sapiens                                                                                                                                                                                                                                                           |
| IgG type                     | Human IgG4 - Kappa                                                                                                                                                                                                                                                     |
| Reactivity                   | Human                                                                                                                                                                                                                                                                  |
| Target                       | LAG3                                                                                                                                                                                                                                                                   |
| Uniprot ID                   | P18627                                                                                                                                                                                                                                                                 |
| Description                  | Anti-LAG3 (relatlimab biosimilar)mAb                                                                                                                                                                                                                                   |
| Delivery                     | In Stock                                                                                                                                                                                                                                                               |
| Storage & Shipping           | Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.                    |
| Background                   | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.                                                                                                                                                                  |
| Usage                        | Research use only                                                                                                                                                                                                                                                      |
| DIMA Disclaimer              | All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement. |





**Figure 1.** ELISA plate pre-coated by 1 µg/ml (100 µl/well) Human LAG3, His tagged protein ([getskuurl sku="PME100512"]) can bind Anti-LAG3 Neutralizing antibody (BME100062) in a linear range of 0.6-100ng/ml.



**Figure 2.** Flow cytometry analysis with Anti-LAG3 mAb 1 µg/ml (relatlimab) on Expi293 cells transfected with human LAG3 (Blue histogram) or Expi293 transfected with irrelevant protein(Red histogram).

